


**PG074120 : BOETTCHER TONSIL SCISSORS**

Curved, Serrated, 7 1/4" (18.5cm)

**PG507453 : ASPIRATING NEEDLE**

Reinforced angular shaft, with protective cover, distal needle 23 gauge (0.6 mm), 14.7 mm long, working length 8 1/4" (21 cm)

**PG507454 : SPARE DISTAL NEEDLE TIP FOR ASPIRATING NEEDLE 507453**

23 gauge (0.6 mm), 14.7 mm long

**PG514430 : BERNAY TRACHEAL RETRACTOR**

Sharp, single prong, 9.5mm depth blade, 6 5/8" (17.0cm)

**PG514431 : BERNAY TRACHEAL RETRACTOR**

Sharp, Double-pronged, 8mm width x 9.5mm depth blade, 6 5/8" (17.0cm)

**PG514433 : BERNAY TRACHEAL RETRACTOR**

Sharp, 4-pronged, 1/2" (13.5mm) width x 9.5mm depth blade, 6 5/8" (17.0cm)



PG514435 : BERNAY TRACHEAL RETRACTOR

Blunt, single prong, 9.5mm depth blade, 6 5/8" (17.0cm)



PG514436 : BERNAY TRACHEAL RETRACTOR

Blunt, Double-pronged, 9.5mm depth blade, 6 5/8" (17.0cm)



PG514437 : BERNAY TRACHEAL RETRACTOR

Blunt, 3-pronged, 9.5mm depth blade, 6 5/8" (17.0cm)



PG514438 : BERNAY TRACHEAL RETRACTOR

Blunt, 4-pronged, 9.5mm depth blade, 6 5/8" (17.0cm)



PO411018 : MOUTH GAG DOYEN 2.5CM ADULT



PO411020 : TONGUE DEPRESSORS DOUBLE STEEL ADULTS



PO411042 : MOUTH GAG JENNING 9CM

PO411043 : MOUTH GAG JENNING 11CM



PO411044 : MOUTH GAG JENNING 13CM



PO411045 : MOUTH GAG JENNING 15CM



PO411055 : COMPLETE MOUTH GAG KILNER DOUGHTY WITH 3 TONGUE DEPRESSORS



PO411207 : TONGUE DEPRESSORS DOUBLE STEEL CHILDREN

